

Comfort Hold™ Thigh & Foot Positioner

The Bionix Comfort Hold Thigh & Foot Positioner is one of the newest products in the line of comfortable patient positioning devices.

Available as either a system or individual units.



With each purchase of this system, Bionix is donating a portion of the proceeds to prostate cancer research.


This system offers multiple indexing options, thus providing one of the most customizable devices in the market.

Features:

- Foam construction for leg comfort.
- Carbon fiber foot board construction for strength.
- Natural positioning of feet, reducing leg fatigue.
- Removable foot cup pads to accommodate shoes-on and shoes-off treatment preferences.
- Attaches to most couch tops.



Entire System.



Sample patient position.



*The products you want,
the service you deserve.*

5154 Enterprise Blvd. - Toledo, OH 43612
phone 800.624.6649 www.BionixRT.com




Comfort Hold Foot Positioner RT-6030-50

- Foam lined foot cups provide additional comfort.
- Multiple indexing locations allowing the foot plate to move anterior or posterior to the patient in the supine position.
- The vertical foot plate will have an option of 5° slightly forward or outward angle at 22.5°.
- Two options for height adjustments and a straight or a splayed angle at 20° for the feet.
- Attaches to most couch tops.



Comfort Hold Thigh Positioner RT-6030-30

- Constructed from a comfortable foam for increased patient comfort.
- Consists of a dual channel riser to support the legs at a 115° bend at the knee.
- 3" riser for additional height.
- Attaches to most couch tops.



*The products you want,
the service you deserve.*

5154 Enterprise Blvd. - Toledo, OH 43612
phone 800.624.6649 www.BionixRT.com

