

# Klarity® Thermoplastic from Bionix

## Klarity Contour and EZ Frames

Klarity Contour and EZ Frames represent the most popular U-Frames style in the industry for several reasons.

- Low Melt Temperature: Reduces set-up time and shrinkage for maximum patient comfort.
- Smooth Edge Comfort: Reduces discomfort around neck line.
- Low Profile Frame: More beam angles can be utilized.
- ID Label
- Fit most head and neck boards in the industry.

### RT-450U Klarity Contour Frames (9" x 10"; 2.5mm)

- Reinforced nose for more support.
- Contoured, perforated patch pattern for easier breathing and increased comfort at the eyes and mouth.

### RT-430UC Klarity EZ Frame (9" x 10"; 3.2mm)

- Thicker for a more rigid mold.
- Rounded edge exclusively from Bionix.

### RT-430W Klarity Profile Frame (10" x 10"; 3.2mm)

- Narrow frame and more thermoplastic.
- Extended version also available.



## Klarity VersaFrames

The VersaFrame and SecureFrames series offers immobilization to the IMRT therapy technique featuring:

- Extended and larger perforations around the eyes, mouth and nose.
- Larger perforations near the head and shoulder to allow for easier breathing.
- NEW Secure frame utilizing an easy-to-use push/pull expansion knob system.
- Fits most industry standard head boards.



Reinforced IMRT SecureFrames are also available.



*The products you want,  
the service you deserve.*

5154 Enterprise Blvd. - Toledo, OH 43612  
phone 800.624.6649 [www.BionixRT.com](http://www.BionixRT.com)



# S-Type Head & Shoulder and Head Only Frames

An alternative to other high-cost S-Type thermoplastic, the Klarity Disposable S-Type frames answer the call.

## Features:

- NEW Easy-to-use push/pull expansion knob system.
- Minimal shrinkage.
- Will not stick to patients or other equipment.
- Low melt temperature.

## 6 Different Styles to choose from:

### Head & Shoulder

- RT-460SB 2.4mm
- RT-461SB 3.2mm
- RT-464SB 3.2mm Reinforced

### Head Only

- RT-460SHB 2.4mm
- RT-461SHB 3.2mm
- RT-464SHB 3.2mm Reinforced



\*US Patent #s: 7,213,597 & D535,024

# Waterbaths for Thermoplastics

## RT-750

Perfect for the compact size of head only thermoplastic. The skillet measures 14.5" x 10.5" x 2.5" deep. 100v only.

## RT-1501 / RT-1502

The model has a spacious design allowing ample room for an 18" x 24" thermoplastic sheet. It also accommodates head and shoulder thermoplastic for easier melting and handling of the product. The hinged and insulated stainless steel lid is designed to open completely or one side only.

## Other Features

- Pre-Set 175° F (79.4°C) on High; 160° F (71.1°C) Low
- Independent heating elements.

Specifications	RT-1501 / RT-1502
Overall Dimensions	28 ½"w x 21"L x 7"H
	72cm x 53cm x 18cm
Interior Tank Dimensions	24 ½"w x 18 ½"L x 5"H
	62cm x 47cm x 13cm
Voltage Requirements	120 Volts
	(220 volts if modified)
Amperage/Hertz	12 amps/60hz
	(with hospital grade plug)
Wattage	720 watts (1440 watt warm-up circuit)



*The products you want,  
the service you deserve.*

5154 Enterprise Blvd. - Toledo, OH 43612  
phone 800.624.6649 www.BionixRT.com

