

# VersaBoard™ Accessories

## Klarity VersaFrames

Reinforced IMRT SecureFrames are also available.



RT-460VHB\*



RT-460VB\*

4 VersaFrame styles  
to choose from:

- RT-460VHB 2.4mm
- RT-461VHB 3.2mm
- RT-460VB 2.4mm
- RT-461VB 3.2mm

2 SecureFrame styles  
to choose from:

- RT-464VHB 3.2mm
- RT-464VB 3.2mm

The VersaFrame mask series are designed specifically for the Bionix VersaBoard. The SecureFrame series feature reinforcements designed with IMRT in mind and fit most industry standard head boards. Both series feature our new T-Pin push/pull locking system.

- **NEW** Secure frame utilizing an easy-to-use push/pull expansion knob system.



*The products you want,  
the service you deserve.*

5154 Enterprise Blvd. - Toledo, OH 43612  
phone 800.624.6649 [www.BionixRT.com](http://www.BionixRT.com)

