

Max3™ Plus Breast Board

CT Friendly and IMRT ready, the **New Carbon Fiber Max3 Plus Breast Board** is 45cm at the widest point. Minimal loose parts and ergonomic design are highlighted by the fully integrated carbon fiber wedge system and secure Velcro closure straps providing ease of use and transport.



The Max3 Breast Board Plus features the following benefits:

- Four integrated elevation carbon fiber wedges at 7°, 12°, 17°, and 23°.
- Indexing Backstop Cushion to prevent patients from involuntarily sliding downward.
- Able to accommodate patients up to 350lbs.
- Compatible with Klarity Head Only VersaFrames.
- Compatible with Silverman head supports.
- Locks down onto most treatment couches.
- Accommodates Butterfly Board, Arm Supports and SecureVac™ Bags
- Large grid panels for increased treatment angles

Breast Board Accessories



Clocking Plate Arm Support
RT-5035-60 Left Support
RT-5035-70 Right Support



RT-5043
Calibrated Arm Cups



RT-460VH
Head Only VersaFrame



*The products you want,
the service you deserve.*

5154 Enterprise Blvd. - Toledo, OH 43612
phone 800.624.6649 www.BionixRT.com



Butterfly Boards



RT-7527 Carbon Fiber Extended Butterfly Board



RT-7526 Carbon Fiber Butterfly Board

All Butterfly Boards include adapter posts to allow for breast board compatibility.
Acrylic Butterfly Boards also available.

SecureFit™ Bars

SecureFit Bars are also available in or order to securely attach the Max3 to your treatment table.



RT-7010 Carbon Fiber SecureFit™ Bar Varian Exact



RT-5010A SecureFit™ Bar Varian Exact



RT-5018 Elekta™ Bar, **RT-5016** Siemens Bar



RT-5012A SecureFit™ Bar II Varian Exact

Indexable sliding bar allows for off-center treatments.

*The products you want,
the service you deserve.*

5154 Enterprise Blvd. - Toledo, OH 43612
phone 800.624.6649 www.BionixRT.com

