

Multi-Chamber SecureVac™ Cushions

The Multi-Chamber SecureVac Cushions are the newest addition to Bionix's SecureVac Cushion Line. The innovative, patented design enables therapists to easily create large full body molds. Each chamber is formed independently, which simplifies the process of full body molds. It's also easy to recreate a section without losing the entire mold. The cushions work well with the Bionix Omni V® SBRT Patient Positioning System, Butterfly Boards or as an independent patient positioning device.

RT-7600 Tri-Vac Cushion*

215 x 70cm, 52L Fill

The Tri-Vac Cushion is designed for full body immobilization.

The cushion has three separate chambers:

- Upper Chamber 65L x 70Wcm, 12L Fill
- Middle Chamber 65L x 70Wcm, 12L Fill
- Lower Chamber 90L x 70cm, 35WL Fill

RT-7625 T-Vac Cushion*

165L x 130cmW, 70L Fill

The T-Vac cushion has two separate chambers and is shaped like the letter "T". The shape allows the SecureVac Cushion to be a popular choice for immobilizing upper extremities. It works extremely well when utilized with a Butterfly Board.

- Upper Chamber 60L x 130cmW, 35L Fill
- Lower Chamber 100L x 70cmW, 30L Fill

RT-7650 Tri-Vac Cushion with Integrated Thigh Bolster*

215 x 70cm, 52L Fill

The Tri-Vac Cushion with Thigh Bolster is the perfect solution for patients that need extra support under their knees. The built-in wedge can easily be positioned anywhere within the cushion's lower chamber.

- Upper Chamber 65L x 70Wcm, 12L Fill
- Middle Chamber 65L x 70Wcm, 12L Fill
- Lower Chamber 90L x 70cm, 35L Fill
- Built-in wedge 6" x 18"L



MR Safe

BUILT-IN ANTIMICROBIAL PROTECTION

*Patented

Multi-Chamber SecureVac™ Cushion Accessories



RT-7101 SecureVac Pump



RT-7527 Carbon Fiber
Extended Butterfly Board



RT-6017 T-Form SecureVac
Box Adapter

*The products you want,
the service you deserve.*

5154 Enterprise Blvd. • Toledo, OH 43612
Phone 800.624.6649 • www.BionixRT.com

