

# IN THIS ISSUE

- ▶ **Bionix Introduces the NEW Prone Breast System**
- ▶ **Order products from the Bionix web site**
- ▶ **Meet Rachel Masters**
- ▶ **Products you want, service you deserve**
- ▶ **Trivia Questions**



Bionix Radiation Therapy  
P.O. Box 935  
Toledo, Ohio 43697-0935

**AFTER ASTRO  
PROMOS INSIDE**

PRST STD  
US POSTAGE  
PAID  
TOLEDO, OHIO  
PERMIT No. 984

## BIONIX® - Products You Want. Service You Deserve.

### New Prone Breast System

Bionix Radiation Therapy unveiled the NEW Bionix Prone Breast System at the 2006 ASTRO Convention. The Prone Breast System is the newest product in the line of comfortable patient positioning devices. The New Prone System features a foam cushion set designed to attach to the Multi-System Base, also used with the Pelvis/Belly System.

- Indexing patient support saddle holds the patient in position for comfort and reproducibility.
- Lightweight design; total weight (including Multi-Treatment Base) is 18 pounds.
- Available complete with Multi-Treatment Base (part # 6250), or as a Cushion Set (part # 6250-10), for those customers with a Multi-Treatment Base from the Pelvis/Belly System.



The NEW Prone Breast System

The new Prone Breast System is a perfect combination of patient comfort and reproducibility. This system is a must have for all facilities currently utilizing prone positioning or that may begin treating breasts in the prone position.


Call today for a demo or to receive a quote on the Bionix Prone Breast System.

#### Key Features & Benefits:

- Utilizes prone positioning for breast treatment to deliver less irradiation to the lung and cardiac tissue.
- Comfortable foam design cradles the patient during treatment cycle.
- Indexing, massage-style face cushion supports the head.
- Reversible, carbon fiber breast bridge holds the contralateral breast comfortably out of the treatment field.

Phone: 1.800.624.6649

### Meet Rachel Masters



Rachel is a recent graduate from the University of Toledo where she completed a double major in marketing and professional sales. Rachel joined the Bionix Radiation Therapy Sales Team in January 2006 as a Customer Service Representative. As a result of her outstanding performance in Customer Service, she is transitioning to Outside Sales where she will serve as the primary contact for a group of Bionix customers.

Fax: 1.800.455.5678

### New Shoulder Suppression

Bionix Radiation Therapy introduces the NEW Premium Shoulder Suppression System. The Premium Shoulder Suppression System features three axes of adjustability, for enhanced patient comfort and precise patient positioning. The Premium Shoulder Suppression System is manufactured from Carbon Fiber to limit attenuation in the treatment fields.

Bionix developed the Premium Shoulder Suppression System for compatibility with the VersaBoard™ and Narrow VersaBoard™.



The NEW Shoulder Suppression System

The VersaBoards offer the latest designs for treating head and neck patients off the couch with IMRT.

**Contact Bionix today for a quote on the new Premium Shoulder Suppression System or for a complete head and neck set-up including the VersaBoard. Call 800.624.6649, email [RTSales@Bionix.com](mailto:RTSales@Bionix.com), or submit a quote request online at [www.BionixRT.com](http://www.BionixRT.com).**

Web: [www.BionixRT.com](http://www.BionixRT.com)

# THE RT REPORT

A NEWSLETTER FROM BIONIX® RADIATION THERAPY

Products You Want. Service You Deserve. BIONIX® RADIATION THERAPY


**Attention! Now Purchase Bionix Products Online!**

Recently Bionix opened the Radiation Therapy Online Store, making it even easier to purchase your consumable products, as you need them. Products available at the Online Store include:

- Klarity® Thermoplastic
- Head & Neck Reusable Frames
- MoldCare® Pillows
- SecureVac™ Cushions & Accessories
- SecureFoam™ two part foam

To maintain a secure website, you will be asked to submit facility and contact information prior to utilizing the Online Store the first time.

Within 24 hours Bionix will send a confirmation of your account and you will be given access to the store. It's that simple. Each return visit to the Online Store simply login with your username and password and order the products you need.



**Products You Want. Service You Deserve.**

Bionix Radiation Therapy is committed to providing the quality products you want with the highest level of service in the industry. The following is a list of the categories of immobilization products Bionix offers along with the service



options that set us apart from the industry.

**Product Lines**

**Head and Neck Immobilization -** Bionix offers a complete line of immobilization equipment for head and neck treatments. This product category features the VersaBoard, designed to treat patients off the couch when treating with IMRT. Accessories for the VersaBoard include the New Premium Shoulder Suppression System with three axes of adjustability.

**Breast and Lung Immobilization -** The Bionix breast and lung product category includes the New Prone Breast System and the Max3 Narrow Breast Board.

**Pelvis and Belly Immobilization -** Bionix features the Pelvis/Belly Plus System for all of your pelvic or abdominal treatments. The Pelvis/Belly Plus System focuses on patient comfort, designed with cushions that cradle the patient during treatment. The Pelvis/Belly Plus utilizes the Multi-Treatment Base, the same base used for the New Prone Breast System.


**SecureVac Cushions -** The Bionix SecureVac Cushions are great for your breast and lung or pelvic set-ups.

**MoldCare Pillows -** MoldCare Pillows, available from Bionix, provide you with the ability to custom mold head cushions for your patient set-ups. MoldCare Pillows index to each of the Bionix head and neck boards for increased reproducibility.



Web: [www.BionixRT.com](http://www.BionixRT.com)

**Klarity Thermoplastic® - Bionix** offers a large variety of Klarity Thermoplastics for use with all your patient positioning needs. Popular options of Klarity Thermoplastics include head only Contour Frames, EZ Frames, Profile Frames and IMRT SecureFrames for head only or head and neck treatment on the VersaBoard.



**Services Offered**  
**Top Level Customer Support -** Bionix Radiation Therapy has developed a service focused sales and customer relations team. The goal of our team is to understand your radiation therapy immobilization needs and match them with quality products. Valuable services offered by Bionix include:


**Product Demo Program -** The Bionix Customer Relations Team will work with your facility to provide you the products needed for proper evaluation. We are confident that once your facility has a chance to use our products, the value they provide will be clear.

Phone: 1.800.624.6649

**In-Service Training Sessions -** The Bionix Sales Team is also available for in-service training sessions at your facility for your entire team. Face to face training allows our representatives to answer any of your questions with the product right in front of you.

**Timely Quoting Process -** Bionix takes pride in providing you with timely price quotes following your request.


Typically, quotes are sent out on the same day as your request, eliminating any delays in your purchase process.



Fax: 1.800.455.5678

**I Know The Answer, Do You?**

Test your knowledge and the knowledge of your colleagues in the Radiation Therapy Department with a five-question quiz. Circulate these questions to each of your colleagues and see who comes out on top.



1. Gala, Winesap, Pippin and Fuji are all varieties of what?
2. What is the correct response to the Irish greeting "Top of the morning to You?"
3. The game of basketball was first played with what type of ball?
4. Who was the first person chosen as the "Times Magazine Man of the Year"?
5. What is the official bird of New Mexico?

The answers to this editions trivia questions can be found by visiting [www.BionixRT.com](http://www.BionixRT.com).

Web: [www.BionixRT.com](http://www.BionixRT.com)

Phone: 1.800.624.6649

Fax: 1.800.455.5678