

T-Form Extremity Immobilizer System

RT-6060 T-Form Extremity Immobilizer System:



Extremity soft tissue sarcomas are some of the most challenging areas to position and immobilize. The T-Form Extremity Immobilizer was designed to provide stability, comfort and high precision, along with the ability to conform to different patient anatomical tumor locations that vary in size and shape.

Features:

- Durable, light-weight, radiolucent T-Form Base Plate
- System specific thermoplastic sheets with reusable handles
- System specific SecureVac Cushions
- SecureFit Bar with Sliding Attachment enables device to be positioned and indexed anywhere on the treatment couch top.
- Locks down onto most couch tops

System Includes:

- 1 - T-Form Base Plate
- 4 - thermoplastic Sheets
- 2 - SecureVac Cushions
- 1 - SecureFit Sliding Attachment
- 1 - SecureFit Lockdown Bar



T-FORM
Radiolucent. Light-weight. Durable.



T-Form Extremity Immobilizer System Accessories



**RT-7136 Extremity
SecureVac Cushion**



**RT-340P Extremity
Thermoplastic Sheets**



**RT-5012A SecureFit™
Bar II Varian Exact**

*The products you want,
the service you deserve.*

5154 Enterprise Blvd. - Toledo, OH 43612
phone 800.624.6649 www.BionixRT.com

