

T-Form Head & Neck Immobilizer

One of the **NEWEST** items available from Bionix!

Bionix has developed a new head immobilization device designed to offer a sensible alternative to carbon fiber head boards.

From imaging and simulation to treatment, the new T-Form Head Immobilizer will be a trusted device in any radiation therapy department.

Features:

- MRI compatible radiolucent material.
- Cost effective alternative to carbon fiber.
- Lighter weight than carbon fiber: Net weight 1lb.
- Fits both disposable and reusable frames.
- Locks down to most treatment couches.

Compatible with these other Bionix devices:

- Carbon Fiber 45° Tilt Attachment
- Prone Attachment
- Lateral Attachment
- Head Only Klarity Frames
- Silverman, Moldcare, and VersaCushion™ Head Supports



RT-7664
45° Tilt Attachment

RT-7029
Prone Attachment

RT-450U
Klarity® Contour
Frame

RT-6200
MoldCare®

RT-7135
VersaCushion™

*The products you want,
the service you deserve.*

5154 Enterprise Blvd. - Toledo, OH 43612
phone 800.624.6649 www.BionixRT.com

