

Butterfly Boards

BREAST & LUNG POSITIONING



RT-7526
Carbon Fiber Butterfly Board



RT-7527
Carbon Fiber
Extended Butterfly Board


Bionix's standard and extended Butterfly Boards are designed to comfortably support patient's arms during treatment procedures of the breast, lung and thorax.

- Available in carbon fiber and acrylic
- Locks down to most couchtops.
- Extended Butterfly Boards offer a T-Grip handle for two axis of adjustability.
- Butterfly Boards easily attach to Bionix Breast Boards.

 Made in the U.S.A.



RT-5526
Acrylic
Butterfly
Board



RT-5527
Acrylic
Extended
Butterfly
Board




RT-5528
Acrylic
Narrow
Extended
Butterfly
Board

Various uses of Butterfly Boards



Butterfly Boards attach to the Bionix Max3 Breast Board.




Flat supine position with T-grip handle and SecureVac Cushion.



Supine position demonstrating frog-leg positioning with SecureVac Cushion.

Accessories



RT-8100 Silverman, Set A-F



RT-7535
SecureVac Cushion
100 x 70, 35L Fill