


Patient Name \_\_\_\_\_

Date \_\_\_\_\_

Area of RX- orientation \_\_\_\_\_

Patient I.D. \_\_\_\_\_



\*For all that index to the base, all readings should be taken from reading the white dot on the side of the apparatus.

### T-Form Treatment Base

Indexed to couch @ \_\_\_\_\_ & \_\_\_\_\_

SecureFit Bar TOP:  A1  B1  C1  D1  E1  F1

SecureFit Bar Bottom:  A2  B2  C2  D2  E2  F2

### Upper Arm Support

\* Index to Base: \_\_\_\_\_ mm

Vertical Height:  A  B  C  D  E  F

Superior to Inferior (0-15 cm) : \_\_\_\_\_ cm

(Left Arm Angle) **Forearm Cups** (Right Arm Angle)

1  2 |  1  2

### Thigh Positioner

\* Index to Base: \_\_\_\_\_ mm

With Riser  Without Riser

### Foot Positioner

Foot Plate Position (0-16 cm): \_\_\_\_\_ cm

(Left Foot) **Cup Position** (Right Foot)

Straight  Straight  
 Splayed  Splayed

Cup Padding:  In  Out

### Thigh Strap (not pictured above)

Shin  Knee  Thigh

Velcro Belt @ \_\_\_\_\_ & \_\_\_\_\_

### Respiratory Belt (not pictured above)

\* Index to Base @ \_\_\_\_\_ mm

Velcro Belt @ \_\_\_\_\_ & \_\_\_\_\_

Pressure Gauge @ \_\_\_\_\_

Insert:  Yes  No

### SecureVac Cushion Used

T-Vac Cushion (T-Shape, 2 Chamber)

Tri-Vac Cushion (3 Chambers)

Other: \_\_\_\_\_

### Turner Arches

#### Short Arch

\* Index to Base @ \_\_\_\_\_ mm

Vertical Index Course Adjust

(S1-S9) \_\_\_\_\_

Respiratory Plate

Fine Adjustment

(10-90) \_\_\_\_\_

#### Tall Arch

\* Index to Base @ \_\_\_\_\_ mm

Vertical Index Course Adjust

(T1-T16) \_\_\_\_\_

Respiratory Plate

Fine Adjustment

(10-90) \_\_\_\_\_