

# All inclusive, sterile, single-use tattoo device

## Accu-Tatt™

The Accu-Tatt is an all inclusive, sterile, single-use tattoo device for patient marking.

Every Accu-Tatt Pouch conveniently contains (1) 1", 18 or (1) 1" 22 gauge safety hypodermic needle and (1) ampoule with black tattoo ink.

The tattoo markings made by Accu-Tatt are both small and permanent. The ergonomic design allows for increased control and assists in accurate patient setups. There are 25 tattoo pouches in every Accu-Tatt box.

### Features:

- Convenience - All inclusive system.
- Safety - Product sterility.
- Single patient use.
- Standardizes the administration of tattoos.
- Increased control.

### Product Offering:

#### **UTRT-3100 Accu-Tatt (18 gauge needle)**

25 pouches/box (Includes 1", 18 gauge safety needle and a black ink ampoule in every pouch)

#### **UTRT-3105 Accu-Tatt (22 gauge needle)**

25 pouches/box (Includes 1", 22 gauge safety needle and a black ink ampoule in every pouch)



MKT-0010-00R01

Phone 800.624.6649 • [www.BionixRT.com](http://www.BionixRT.com)  
5154 Enterprise Blvd., Toledo, Ohio 43612

