

# Comfort Hold™ Thigh & Foot Positioner

The Comfort Hold Thigh and Foot Positioner is a customizable device offering multiple indexing options for the lower extremities. This versatile product allows natural positioning of the feet and legs, reducing leg fatigue.

The Comfort Hold Thigh and Foot Positioner is available as a system or can be purchased as individual units.

## PBRT-6030 Comfort Hold Thigh & Foot Positioner



Bionix will donate a portion of the proceeds from the purchase of this system to prostate cancer research.

### Features:

- Foam construction wedges
- Lightweight radiolucent T-Form Base
- Natural positioning of feet, reducing leg fatigue
- Attaches to most couch tops
- Removable foot cup pads



## Comfort Hold™ Thigh & Foot Positioner Accessories



**SVRT-7535**  
SecureVac Cushion



**SVRT-7625**  
T-Vac Cushion

Phone 800.624.6649 • [www.BionixRT.com](http://www.BionixRT.com)  
5154 Enterprise Blvd., Toledo, Ohio 43612



## Comfort Hold Foot Positioner

### PBRT-6030-30-03

- Multiple indexing locations allows the light, radiolucent foot plate to move superior to inferior
- Offers two height adjustments options and a straight or 20° splayed angle for the feet
- Removable foot cup pads to accommodate shoes-on and shoes-off treatment preferences
- Locks down onto most treatment couches




## Comfort Hold Thigh Positioner

### RT-6030-50-02

- Constructed from durable foam to aid in patient comfort
- Dual channel riser to support the legs at a 115° bend at the knee
- 3" removable riser for multiple height options
- Locks down onto most treatment couches

Shown with  
3" riser.



Without 3" riser.

