

The Bionix Radiation Therapy *T-Form Head and Neck Immobilizer*


The T-Form Head & Neck Immobilizer is made out of a composite material that is radiolucent, light-weight and durable. The T-Form Head & Neck Immobilizer fits most U-Frame thermoplastic masks and interfaces with various attachments to provide positioning options.

Features:

- Carbon Fiber-like attenuation
- Capable of supine, prone or lateral treatment positions
- Fits most U-Frame thermoplastic masks

Compatible with these other Bionix devices:

- Carbon Fiber 45° Tilt Attachment
- Prone Attachment
- Lateral Attachment
- Head Only Embrace Frames
- Silverman, MoldCare, and VersaCushion™ Head Supports



© 2018 Bionix Radiation Therapy • MKT-0048-00 • R02

