


# The Bionix Radiation Therapy *Dry Oven*

The BRT Dry Oven heats thermoplastic material using safe heat technology, while reducing the risk of cross-contamination. An intuitive digital display allows the dry oven to heat the thermoplastic material at an optimal temperature for molding, ensuring maximum patient comfort during application.



## Features:

- Compact, space-saving design, ideal for crowded departments
- Eliminates the inconveniences of using water
- Non-stick surface
- 3" Touch Display for Temperature adjustments (Celsius or Fahrenheit)
- Insulated case
- Built-in timer function
- Outside Dimensions: 79cm x 58cm x 13cm
- Inside Dimensions: 66cm x 51.5cm x 5.8cm

## Product Offering:

**KLSP-1601-D - Dry Oven**

MKT-0018-00R01

Phone 800.624.6649 • [www.BionixRT.com](http://www.BionixRT.com)  
5154 Enterprise Blvd., Toledo, Ohio 43612

